



# Designing selective irreversible inhibitors



# Designing selective irreversible inhibitors

Irreversible inhibition involves an initial binding step driven by affinity and a time-dependent inactivation step driven by covalent bond formation. As  $IC_{50}$  does not reflect both steps, its use may fail to define selectivity.

Kinetic characterization allows for the determination of the rate ( $k_{inact}$ ) and efficiency ( $k_{inact}/K_I$ ) of inactivation as well as the potency ( $K_I$ ).

This data is more useful for the design of selective covalent inhibitors than  $IC_{50}$  values, since it allows for the improvement of both inhibitor binding and covalent bond formation. Modification of any part of an inhibitor may improve the rate of inactivation ( $k_{inact}$ ) at the expense of inhibitor binding ( $K_I$ ) or vice versa.

In order to minimize potential off-target effects, the right balance between efficiency of inactivation and selectivity is needed.

A recent study demonstrated how  $IC_{50}$  inadequately reflects the selectivity profile of the BTK irreversible inhibitors ibrutinib and acalabrutinib. Although acalabrutinib seems to be more selective than ibrutinib when comparing  $IC_{50}$  alone<sup>3</sup>, the assessment of selectivity using  $K_{inact}/K_I$  ratios corroborate clinical data, demonstrating similar safety profiles between the therapies.<sup>1-3</sup>

## COVALfinder<sup>®</sup> helps you to modulate the therapeutic window and safety profile of irreversible drugs

### Features and benefits

- Accurate.
- Robust and reproducible.
- Sensitive.
- Broad dynamic range.
- Rapid turnaround.
- Standard known inhibitor in each plate.

### Applications

- Modify compound binding ( $K_I$ ) and reactivity ( $k_{inact}$ ).
- Modulate the therapeutic index and safety profile of irreversible drugs.
- Distinguish between one- or two-step irreversible mechanisms.
- Understand PK/PD disconnects.

## How COVALfinder<sup>®</sup> assays work

COVALfinder<sup>®</sup> is a highly tuned TR-FRET Kinetic Assay. Real-time binding of an active-site directed fluorescent probe is detected using a labeled anti-tag antibody, which binds to the target of interest. The binding of the probe to the target increases the TR-FRET signal, whereas the displacement of the probe with a compound decreases the TR-FRET signal.

### Assay process

#### 1. Performance:

- Pre-incubation of the target of interest, a fluorescent probe and labelled antibody.
- 20-point 1.75-fold serial dilutions in duplicate of test compounds and a reference compound.

- Up to 48 total binding and non-specific binding controls.
2. **Detection:** Reaction is monitored over time at room temperature.
  3. **Analysis:**
    - Fitting of specific TR-FRET signals to a single-exponential equation to obtain  $k_{obs}$ , the first-order rate constant for the interconversion between the initial and final degree of binding.
    - Production of dose-response curves at each time point in order to generate  $IC_{50}$  values, allowing easy inspection of time dependency.
    - Determination of  $K_I$ ,  $k_{inact}$ ,  $k_{inact}/K_I$  values and the mechanism of inhibition (one- or two-step) using a secondary plot of  $k_{obs}$  versus compound concentration and a statistical analysis.
  4. **Quality control:** S/B, Z-value and MSR (within and between-run variability) for each assay.

## Selectivity profiling of irreversible EGFR inhibitors using COVALfinder<sup>®</sup>

EGFR is a receptor tyrosine kinase, member of the ErbB family together with HER2, HER3 and HER4. EGFR overexpression is frequent in many types of human malignancies, including non-small cell lung cancer (NSCLC). Afatinib, Dacomitinib and Neratinib are quinazoline based irreversible pan-ErbB drugs. Afatinib and Dacomitinib have been approved for the treatment of EGFR mutated NSCLC and Neratinib for the treatment of HER2-Positive breast cancer.

Here we study the selectivity profile of three irreversible EGFR inhibitors using COVALfinder<sup>®</sup>. The binding profiles for all the irreversible inhibitors showed a slow onset of

EGFR, HER2 and HER4 inhibition (representative data for HER2 inactivation by Dacomitinib is shown in Fig.1, A). Each inhibitor with the three kinases tested also




Figure 1. **Characterization of the kinetics of HER2 inactivation by the irreversible inhibitor Dacomitinib.** (A) Progress curve of HER2 incubated with increasing compound concentrations. (B) Dependence of  $k_{obs}$  on Dacomitinib concentration. (C) Dose-response curves over time. (D)  $IC_{50}$  values over time.

displayed the characteristic time-dependent reduction in  $IC_{50}$  due to the formation of the covalent complex over time (Fig.1, C and D).

To confirm that these inhibitors were indeed binding irreversibly to EGFR, HER2 and HER4, the dependence of  $k_{obs}$  on inhibitor concentration was evaluated (Fig.1, B). A two-step, time-dependent inhibition was observed. This behavior is consistent with either slow onset reversible inhibition (the slower phase representing the rate of target isomerization) or irreversible inhibition (the slower phase representing target inactivation). However, in all cases the y-intercept was zero, which is only consistent with irreversible binding. The  $k_{inact}$ ,  $K_I$ ,  $k_{inact}/K_I$  and  $T_{1/2}^{\infty}$  (half-life for inactivation at infinite concentration of inhibitor) values obtained from this analysis are presented in Table 1.

## Inactivation kinetics profiling is crucial for selectivity assessment

As can be seen in Table 1, if we take into account the  $IC_{50}$  values to calculate the selectivity ratio, all the inhibitors are pan-ErbB drugs. However, using  $k_{inact}/K_I$  the selectivity profile of Afatinib and Dacomitinib dramatically changes. Afatinib and Dacomitinib binds preferentially to EGFR and HER4 whereas Neratinib is a pan-ErbB inhibitor.




Figure 2. **Covalency quadrant scale.** Classification of covalent binders by target comparing the rate of covalent bond formation ( $k_{inact}$ ) against affinity ( $K_I$ ).

The covalency quadrant (Fig.2) clearly shows that the selectivity of Afatinib and Dacomitinib towards HER2 is due to a lower affinity ( $K_I$ ) and reactivity ( $k_{inact}$ ). These results may explain the efficiency of Afatinib and Dacomitinib in the treatment of EGFR mutated NSCLC and the efficiency of Neratinib in HER2-Positive breast tumors.

Compound	Target	$K_I$ (nM)	$k_{inact}/K_I$ ( $M^{-1}s^{-1}$ )	$k_{inact}$ ( $s^{-1}$ )	$T_{1/2}^{\infty}$ (min)	$IC_{50}$ (nM)	$IC_{50}$ selectivity	Efficiency selectivity
Dacomitinib	EGFR	0.2	$1.6 \times 10^7$	$3.8 \times 10^{-3}$	3.1	6.0	1	1
Afatinib		0.1	$4.6 \times 10^7$	$6.0 \times 10^{-3}$	1.9	0.5	1	1
Neratinib		1.2	$4.4 \times 10^6$	$5.3 \times 10^{-3}$	2.2	92	1	1
Dacomitinib	HER2	8.7	$1.1 \times 10^5$	$9.3 \times 10^{-4}$	12.4	46	8	150
Afatinib		5.5	$5.0 \times 10^5$	$2.8 \times 10^{-3}$	4.2	14	28	91
Neratinib		4.4	$1.2 \times 10^6$	$5.3 \times 10^{-3}$	2.2	59	1	4
Dacomitinib	HER4	7.5	$8.4 \times 10^5$	$6.3 \times 10^{-3}$	1.8	74	12	19
Afatinib		6.3	$2.0 \times 10^6$	$1.3 \times 10^{-2}$	0.9	1.0	2	23
Neratinib		1.5	$7.6 \times 10^6$	$1.1 \times 10^{-2}$	1.0	NA		0.6

Table 1. **Kinetic analysis of irreversible inhibition of ErbB family.** Calculation of the selectivity ratio is based on the  $IC_{50}$  values reported in the literature and  $k_{inact}/K_I$  values obtained with COVALfinder® for off-targets relative to EGFR. NA: no data available in the literature. Empty fields indicate that the selectivity ratio cannot be calculated.

## References

1. Ayah A. et al. (2020) Advances in covalent kinase inhibitors. *Chem Soc Rev.* 49(9):2617-2687.
2. Strelow J.M. et al. (2017) A Perspective on the Kinetics of Covalent and Irreversible Inhibition. *SLAS Discov.*
3. Hopper M. et al. (2020) Relative Selectivity of Covalent Inhibitors Requires Assessment of Inactivation Kinetics and Cellular Occupancy: A Case Study of Ibrutinib and Acalabrutinib. *J Pharmacol Exp Ther.* 372(3):331-338.
4. Engelman J.A. et al. (2007) PF00299804, an irreversible pan-ERBB inhibitor, is effective in lung cancer models with EGFR and ERBB2 mutations that are resistant to gefitinib. *Cancer Res.* 67(24):11924-32.
5. Rabindran S.K. et al. (2004) Antitumor activity of HKI-272, an orally active, irreversible inhibitor of the HER-2 tyrosine kinase. *Cancer Res.* 64(11):3958-65.
6. Li D. et al. (2008) BIBW2992, an irreversible EGFR/HER2 inhibitor highly effective in preclinical lung cancer models. *Oncogene.* 27(34):4702-11.