

Evaluation of covalent BTK inhibitors using COVALfinder[®]

Ana Corrionero, Niall Prendiville, Jose A. Rodriguez, Patricia Alfonso.
Enzymlogi, Qube Technology Park, Madrid, Spain

De-risking covalent BTK drug discovery

BTK is a prominent therapeutic target for hematologic cancers and an attractive target for treating autoimmune diseases. BTK inhibitors can be divided into 2 types: covalent irreversible and non-covalent reversible inhibitors.

Currently, 6 approved BTK inhibitors (Ibrutinib, Acalabrutinib, Zanubrutinib, Tirabrutinib and Orelabrutinib) target the kinase domain of BTK, forming a covalent bond with Cys481. Remibrutinib is another irreversible inhibitor that exhibits a good kinase selectivity due to binding to the non-phosphorylated form of BTK and is being evaluated in the clinic for urticaria and asthma. Recently, a hybrid inhibitor with the ability to establish a reversible covalent bond with Cys481 and temporarily inactivate BTK (Rilzabrutinib) has entered phase 3 clinical trials for the treatment of pemphigus and immune thrombocytopenic purpura^{1,2}.

Here, we present the kinetic characterization of covalent BTK inhibitors using COVALfinder[®] as a successful approach to identify novel covalent drugs, interpret SAR and modulate the safety profile of irreversible drugs.

How COVALfinder[®] assays work

COVALfinder[®] is a robust kinetic platform that delivers all key inactivation (k_{inact} , k_{inact}/K_I , $T_{1/2}^{CO}$) and affinity (K_I) parameters for irreversible binders.




Fig 1. COVALfinder[®] method.

FEATURES

- k_{inact} , K_I , k_{inact}/K_I and $T_{1/2}^{CO}$.
- Robust, reproducible and precise.
- Broad dynamic range.
- Activated and non-activated targets.
- Rapid turnaround time.

APPLICATIONS

- Discrimination between irreversible and reversible drugs.
- Selectivity profiling for irreversible drugs.
- Understanding PK/PD disconnects.

Differentiation of conformation-specific BTK irreversible drugs by inactivation kinetics

The crystal structure of Ibrutinib, Acalabrutinib, Zanubrutinib and Remibrutinib shows the stabilization of different inactive states of BTK^{3,4}. Ibrutinib, Acalabrutinib and Zanubrutinib occupy the back pocket of BTK kinase domain when the C-helix is rotated outwards (C-helix out) and the activation loop adopts the DFG-in conformation. In contrast, Remibrutinib occupies the H3 pocket, which is present when the activation loop Tyr551 is unphosphorylated. Consequently, Remibrutinib blocks the activation loop phosphorylation and stabilize it in the DFG-out conformation⁵.

- Ibrutinib, Acalabrutinib and Zanubrutinib inactivate similarly the phosphorylated and non-phosphorylated forms of BTK, indicating that the inhibitors may occupy the back pocket present in both forms.
- Whereas the inactivation efficiency (k_{inact}/K_I) of Remibrutinib is 140-fold higher with the non-phosphorylated form. This may explain the exquisite selectivity profile of this inhibitor⁵.




Figure 2. Characterization of the kinetics of non-phosphorylated BTK inactivation by the covalent-irreversible inhibitors Ibrutinib, Zanubrutinib, Acalabrutinib and Remibrutinib using COVALfinder[®]. (A) Progress curve of BTK incubated with increasing compound concentrations. (B) Dependence of k_{obs} on inhibitor concentration.

Compound	Phosphorylated BTK				Non-phosphorylated BTK			
	K_I (nM)	k_{inact}/K_I ($M^{-1}s^{-1}$)	k_{inact} (s^{-1})	$T_{1/2}^{CO}$ (min)	K_I (nM)	k_{inact}/K_I ($M^{-1}s^{-1}$)	k_{inact} (s^{-1})	$T_{1/2}^{CO}$ (min)
Ibrutinib	13	1.1×10^6	1.5×10^{-2}	0.8	14	7.2×10^5	1.0×10^{-2}	1.1
Zanubrutinib	37	4.6×10^5	1.7×10^{-2}	0.7	55	2.7×10^5	1.5×10^{-2}	0.8
Acalabrutinib	138	5.6×10^4	7.7×10^{-3}	1.5	118	5.8×10^4	6.8×10^{-3}	1.7
Remibrutinib		1.6×10^4			1.6	2.3×10^6	3.6×10^{-3}	3.2

Table 1. Kinetic analysis of the irreversible inhibition of BTK forms. K_I , k_{inact} , k_{inact}/K_I and $T_{1/2}^{CO}$ of Ibrutinib, Zanubrutinib, Acalabrutinib and Remibrutinib for the phosphorylated and non-phosphorylated form of BTK are shown.

Ibrutinib and Acalabrutinib selectivity profile

Irreversible inhibition involves an initial binding step driven by affinity and a time-dependent inactivation step driven by covalent bond formation. As IC_{50} does not reflect both steps, its use may fail to define selectivity. In order to minimize potential off-target effects, the right balance between efficiency of inactivation and selectivity is needed.

- Ibrutinib irreversibly inhibits kinases which carry a Cys at the same position of BTK: ITK, EGFR, HER2, HER4.
- Acalabrutinib preferentially inactivates HER2 and HER4. Selectivity towards EGFR is due to a lower affinity (K_I) and reactivity (k_{inact}).




Figure 3. A) Covalency quadrant scale. Classification of Ibrutinib and Acalabrutinib by target comparing the rate of covalent bond formation (k_{inact}) against affinity (K_I). B) Kinetic analysis of irreversible inhibition of BTK, HER4, HER2, EGFR and ITK by Ibrutinib and Acalabrutinib. Calculation of the selectivity ratio is based on each parameter for off-targets relative to BTK. Ratios >10 are highlighted in yellow, while those >100 are highlighted in green. Red indicates a similar profile.

In addition, Ibrutinib but not Acalabrutinib, inhibits reversibly several other kinases.

- In all cases, off-target occupancy and dissociation is rapid (seconds to few minutes).
- Ibrutinib maintain activity against BTK owing to prolonged covalent engagement, while the transient reversible inhibition of off-targets is minimized.




Figure 4. A) Kinetic analysis of the reversible inhibition of Ibrutinib with all kinases that have shown a K_d value lower than 100 nM using the KINETICfinder[®] platform. Residence time values within the range of seconds are highlighted in green, while those within the range of minutes are highlighted in yellow. B) Kinetic curves of Ibrutinib for two representative kinases obtained with KINETICfinder[®].

Analysis of BTK reversible-covalent drugs

COVALfinder[®] and KINETICfinder[®] platforms has been used to provide an in-depth understanding of the binding mechanism of the covalent-reversible inhibitor Rilzabrutinib.

- Our results confirm that BTK-Rilzabrutinib complex formation takes place in reversible sequential events.
- Rilzabrutinib binding is the initial event and the covalent bond formation with Cys481 of BTK the second.
- The high affinity of Rilzabrutinib (K_d^* : 0.14 nM) is a result of the extraordinary stability of the BTK-Rilzabrutinib complex (k_{off}^* : $2.0 \times 10^{-5} s^{-1}$, residence time: 812 min).





Figure 5. Characterization of the kinetics of BTK inactivation by the covalent-reversible inhibitor Rilzabrutinib using COVALfinder[®]. (A) Progress curve of BTK incubated with increasing compound concentrations. (B) Dependence of k_{obs} on Rilzabrutinib concentration. (C) Dose-response curves over time. (D) IC_{50} values over time.

Discussion & Significance

- The analysis of K_I , k_{inact} and k_{inact}/K_I of irreversible BTK inhibitors using COVALfinder[®] reveals that Ibrutinib, Acalabrutinib and Zanubrutinib inactivate similarly the phosphorylated and non-phosphorylated forms of BTK. These results show that the inhibitors may occupy the back pocket present in both forms.
- Our results confirm that Remibrutinib occupies the H3 pocket, which is only present in BTK when the activation loop Tyr551 is unphosphorylated: the k_{inact}/K_I obtained with the COVALfinder[®] assay of non-phosphorylated BTK is 140-fold higher than the phosphorylated BTK.
- The evaluation of k_{inact}/K_I of Ibrutinib and Acalabrutinib towards several kinases which carry a Cys at the same position of BTK displays that Acalabrutinib is more selective.
- Our kinetic platform for reversible binders KINETICfinder[®] shows that Ibrutinib transiently inhibits several other kinases, suggesting that the risk of undesired side effects would be minimized.
- COVALfinder[®] enables a detailed characterization of both covalent-reversible and irreversible drugs.
- Rilzabrutinib is a covalent-reversible inhibitor that binds to BTK in two-reversible steps with a fast association rate ($3.3 \times 10^5 M^{-1}s^{-1}$) and very long residence time (812 min).